


# How Affinity Translates into a Selective Advantage: A Comparison of Phage and Yeast Display for Rabbit Mass Humanization

Michael Tomaszowski, Peter Milovnik, Jacob Glanville (Distributed Bio), Remko Griep, Ahmad Trad and Vera Molkenhain

## Introduction

AbCheck's Mass Humanization provides highly diverse human antibodies from immunized non-transgenic rabbits in shorter time compared to hybridoma screens, followed by CDR grafting and antibody engineering. Mass Humanization libraries are generated by batch cloning of the *in vivo* generated CDR3 diversity into human frameworks with designed CDR1/2 diversity, and human antibodies are selected from these libraries. Here, we directly compared Phage panning under different conditions and FACS sorting of a Yeast Display library with respect to their efficiencies in the selection of human high affinity binders from such libraries. The analysis of sequences obtained from both methods revealed overlapping findings, but also significant differences regarding the overall sequence diversity and binding affinity.

**Figure 1. Rabbit Mass Humanization (RMH) libraries**




- 100% human VH/VL framework regions
- No biophysical liabilities
- Designed diversity in CDR1/2
- Combined with CDR-H3 and CDR-L3 diversity from immunized rabbits

**Figure 2. High hit rate and large sequence diversity obtained from low stringency Phage selections**

Results of Sanger sequencing of 300 hits:


- >100 different CDR-H3 sequences
- 21 CDR-H3 sequence families
- Multiple CDR-L3 sequences per CDR-H3
- Varying CDR1/2 mutations

Two enriched CDR-H3 sequence families: Consensus sequence of the two enriched CDR-H3 sequence families:




**Figure 3. Off-rate Phage selection combined with NGS reveals the selective advantage of CDR-H3 sequences**

Relative frequencies of CDR-H3 sequences after three phage panning rounds under varying conditions:




**Figure 4. FACS sorting of a Yeast Display library yields a 100% binding population in 4 sorting rounds**



**Figure 5. FACS sorting at low antigen concentration reveals selective advantage**

Distribution of CDR-H3 sequences found in Sanger sequencing of hits after FACS sorting:



**Figure 6. Properties of selected sequences:**

Sequence	IgG KD	% of V-domain identical to human germline
RMH-yeast 1	5,45E-11 M	93
RMH-yeast 2	1,60E-10 M	93
RMH-phage	1,49E-10 M	92

## Conclusion:

- Phage panning and Yeast FACS sorting are both well suited to identify high affinity (sub-nanomolar) binders from a displayed Rabbit Mass Humanization library
- Based on sequencing, Phage panning yielded a broader sequence diversity which is likely to lead to an extended epitope spectrum and affinity panel of hit sequences compared to Yeast FACS sorting
- One particular CDR-H3 sequence group was selectively enriched in both Phage and Yeast with a significantly better enrichment of high affinity binders out of a highly diverse humanized immune repertoire after Yeast FACS sorting